thermo scientific

# Nicolet iZ10

FTIR instrument accessory

with OMNIC Paradigm software

## SITE & SAFETY GUIDE

269-3648 00

Revision A


March 2024

**ThermoFisher**  
SCIENTIFIC

---

# Set up your new instrument

[https://knowledge1.thermofisher.com/Molecular\\_Spectroscopy/Instrument\\_Setup](https://knowledge1.thermofisher.com/Molecular_Spectroscopy/Instrument_Setup)

---

# Nicolet iZ10 Site & Safety Guide

## MANUAL CONVENTIONS

Safety precautions and other important information use the following format:

### DANGER

**Avoid hazard.** Indicates a hazardous situation which, if not avoided, will result in serious injury or death.

### WARNING

**Avoid hazard.** Indicates a hazardous situation which, if not avoided, could result in serious injury or death.

### CAUTION

**Avoid hazard.** Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

### NOTICE

Follow instructions with this label to avoid damaging the system hardware or losing data.


**Note** Contains helpful supplementary information.

**Tip** Provides information that can make a task easier.

# Site preparation and intended use

The following table lists some of the safety symbols and their indications that may appear in the user documentation.

Symbols	Indication
	This is a mandatory action symbol. It is used to indicate that an action shall be taken to avoid a hazard.
	This is a prohibition symbol. It is used to alert the user to actions that shall not be taken or shall be stopped.
	This is the general warning sign. Failure to heed the safety precautions could result in personal injury. When you see this symbol on your instrument, consult the user documentation for safe operation.
	 This is a sign to avoid shock hazard. If you see either of these symbols, there is a risk of electrical shock in the vicinity. Only qualified persons shall perform the related procedures.
	<b>Avoid fire hazard.</b> Do not test flammable or explosive samples. Read and follow the associated instructions carefully.

Symbols	Indication
	<p><b>Avoid toxic inhalation.</b> Materials such as hydrochloric acid, hydrofluoric acid, and phosgene are highly toxic. If you plan to regularly use solvents containing halogenated hydrocarbons, be sure your work area is properly ventilated.</p>
 	<p><b>Avoid burn injury.</b> These symbols alert you to hot surfaces. Read and follow the associated instructions carefully.</p>
	<p><b>Avoid explosion hazard.</b> These symbols alert you to the risk of a possible explosion. Never use a flammable, combustible, or toxic gas to purge this instrument. Heat from the source, or from laser absorption, may ignite flammable gases or reactive materials in purge gas.</p>
 	<p><b>Avoid eye injury.</b> If you see these symbols, there is a risk of exposure to laser radiation in the vicinity. Read and follow the associated instructions carefully.</p>
	<p><b>Avoid pinch hazard.</b> If you see this symbol, there is a risk of pinch, crush, or impact hazards in the vicinity.</p>



Symbols	Indication
 	<p><b>Avoid risk of personal injury.</b> If you see either of these symbols, use proper lifting techniques when moving the spectrometer or other system components.</p>
   	<p><b>Avoid freeze burns.</b> If you see any of these symbols, there is a risk of skin burns or eye damage. Protective equipment required. Follow standard laboratory safety practices.</p>
	<p><b>Avoid chemical burns.</b> This symbol alerts you to possible skin irritation. Wear gloves when handling toxic, carcinogenic, mutagenic, corrosive, or irritant chemicals. Use approved containers and proper procedures to dispose of waste.</p>

Table 2-1: Electrical symbols

Symbol	Description
~	Alternating current
⏏	Earth terminal or ground
—	Direct current

Symbol	Description
	Protective conductor terminal
	Frame or chassis terminal
	Fuse
	Power on
	Power off

## CAUTION

**Avoid personal injury.** If this equipment is used in a manner not specified in the accompanying documentation, the protection provided by the equipment may be impaired.

**Avoid personal injury.** Perform only those procedures described in the documentation. If there are other problems, contact us. Any other service must be performed by trained personnel.

**Avoid personal injury.** Do not remove the cover of the instrument. All service to the instrument must be performed by trained personnel.

**For Research Use Only. This instrument or accessory is not a medical device and is not intended to be used for the prevention, diagnosis, treatment, or cure of disease.**

## Operating your accessory safely

This safety guide contains a summary of the safety precautions that must be followed when using a Thermo Scientific Nicolet™ iZ10 accessory. Using the accessory in a manner not specified in any of the documentation that came with the instrument may create a hazardous condition. Each person who will be using this instrument should read this manual. Safety information is also included in your spectrometer, accessory, computer, and software documentation.

---

In many cases, safety information is displayed on the instrument itself. The safety label illustrations in the documentation that came with your instrument show the locations of the safety labels. Should any of these labels become loose or unreadable, use the information in the front of this manual to contact us to obtain new ones.

### INGRESS PROTECTION RATING

This product is rated IP20.

### WHEN THE INSTRUMENT ARRIVES

- Check the exterior of the shipping box for signs of damage. If damage is apparent, contact us or your local distributor for instructions.
- Move the shipping box to the installation location at least 24 hours before installation.
- Keep the instrument upright while moving it to the installation location.

### FIRE SAFETY AND BURN HAZARDS

The following symbols appear inside the instrument to warn you about hot surfaces in the vicinity of the symbol.

To avoid a burn injury and the risk of fire or explosion, follow these guidelines:

- Do not test flammable or explosive samples.
- Use only nitrogen or dried air to purge your instrument.
- If there is an infrared source in the instrument, do not touch the source housing; it may be very hot.
- Do not place anything on the electronics cover.
- After you turn off your instrument, wait 15 minutes before you replace components.
- Never block any of the vents on an instrument or its power supply.
- Use exact replacements for power supplies.

The following precautions provide additional information about fire safety and preventing burn hazards.

#### CAUTION

**Avoid hazard.** If there is an infrared source in the instrument, the source housing may be very hot. It stays hot for up to 15 minutes after you turn off your instrument. To avoid being burned or starting a fire, do not touch the housing with your body or any flammable object until it has cooled. The documentation that came with your system shows the location of the infrared source inside your instrument.

## DANGER

**Avoid hazard.** Never use a flammable gas to purge this instrument. The purge gas must be free of oil and other reactive materials. Heat from the source or from laser absorption may ignite flammable gases or reactive materials in purge gas.

## CAUTION

**Avoid hazard.** Do not place anything on top of the vents at the back of the instrument. Electronic components can overheat if these vents are blocked.

### MOVING PARTS PRECAUTION

Some instruments contain moving parts. If moving parts present a hazard, a label alerts you to the hazard.

### TEMPERATURE AND HUMIDITY

- The instrument is designed for indoor use in normal operating environments (Pollution Degree 2)
- Your instrument operates reliably at temperatures between 16 °C and 27 °C (60 °F and 80 °F). It may briefly be exposed to temperatures as low as 10 °C (for example, during transport) without degradation of its safety. Temperature changes may result in drift in the system response.
- Maintain humidity in the range of 20% to 80% non-condensing. Your instrument must be protected from excessive humidity, since the beamsplitters, detectors, and other components may corrode with exposure to moisture. To accomplish this, some instruments are sealed and desiccated; other instruments require purging.
- When stored unpackaged or when operating the instrument in ambient conditions above 40% RH, it is highly recommended that the instrument be purged with dry air or nitrogen to protect the internal components.
- Operation or storage with desiccant only in a high humidity environment may require desiccant replacement every 14 days depending on the conditions.
- Do not rely on the humidity indicator to determine when the desiccant has become saturated; use the System Status table in OMNIC Paradigm software instead.
- Purge the system with clean dry air or nitrogen.
- Avoid rapid changes in temperature that may cause condensation.

---

## ALTITUDE

Your instrument is rated for altitudes of up to 2,000 m (6,500 ft).

## STORAGE

When stored in the original shipping container, the instrument can be exposed to temperatures from -20 °C to 60°C without damage to the instrument. Maximum humidity for storage is 85% RH, non-condensing. The storage time for the instrument at high humidity in its original packaging ranges between 60 and 70 days. Storage at lower ambient humidity will extend this time.

## VIBRATION

- The instrument will perform better in a mechanically stable environment.
- Keep instrument away from machinery that may vibrate the floor.
- Minimize or eliminate acoustic noise and vibration wherever possible.
- Consider placing instrument on a marble top table or counter.

Floor vibration or acoustical noise from heavy manufacturing equipment, computer equipment, or other sources will not damage the system, but it can affect performance and spectral quality.

## MAGNETIC AND ELECTRIC FIELDS

- Place instrument at least 5.5 m (18 ft) away from magnetic fields.
- Minimize or eliminate exposure to magnetic fields wherever possible.
- Some wireless devices may also affect instrument performance. If this type of interference is suspected, move all wireless devices at least 2.0 m (6.5 ft) away from the instrument.

## USING LIQUID NITROGEN

Some detectors must be cooled with liquid nitrogen before use. The following symbols remind you to wear protective clothing when using liquid nitrogen.

## WARNING

**Avoid hazard.** Liquid nitrogen is extremely cold and therefore potentially hazardous. When filling the detector dewar, be careful not to allow the liquid nitrogen to contact your skin. Wear protective gloves and splash-proof goggles and follow standard laboratory safety practices.

To avoid hazardous equipment damage or contact with liquid nitrogen, make sure any dewar or container you use to hold liquid nitrogen can do so safely without breaking.

## LIFTING OR MOVING THE INSTRUMENT

Your instrument weighs 23 kg (50 lb) or more. To avoid risk of injury, use proper lifting techniques when lifting or moving the instrument or other system components.

## CAUTION

**Avoid hazard.** Your instrument cannot be lifted safely by a single person. Lifting or moving this instrument requires two people. Be sure to grasp with both hands under the base on both sides of the instrument.

## SPECIFICATIONS AND WORKSPACE CONSIDERATIONS

**Instrument weight:** 23.0 kg (50.7 lb)

**Dimensions:** 57.15 cm (L), 54.6 cm (W), 31.75 cm (H)

If the optional sample compartment expansion kit is installed, the length (L) is increased by approximately 5.0 cm.

**Computer:** Plan for location of the computer, monitor and keyboard. A standard USB A-B data cable is required to connect the instrument to the computer.

**Clearance:** Leave ample space near and behind the system to allow access to cable and purge connections. The presence of an additional space next to the instrument makes it easier to access, speeding up and facilitating maintenance and operation. Plan for the location of the spectrometer as well, and ensure that the workstation is long enough to accommodate both instruments.

---

# Electrical safety

Your instrument was designed with protective covers to prevent exposure to dangerous voltage and other electrical hazards. If you see either of the following symbols on your instrument, there is a risk of electric shock in the vicinity of the symbol:

Figure 2-1: Shock hazard symbols

## DANGER

**Avoid shock hazard.** To avoid injury, only a qualified person using the appropriate measuring device should check the line voltage, current, or frequency.

## WARNING

**Avoid hazard.** Only our trained and certified service representatives should attempt to service a component that carries either of these symbols.

## WARNING

**Avoid hazard.** If a protective cover on a system component appears damaged, turn off the system and secure it against any unintended operation. Always examine the protective cover for transport stresses after shipping.

## WARNING

**Avoid shock hazard.** Even after this instrument has been disconnected from all voltage sources, capacitors may remain charged for up to 30 seconds and can cause an electrical shock.

## WARNING

**Avoid shock hazard.** Do not allow liquid to run into the power supply. This will result in electric shock.

## NOTICE

Sample compartment windows can be damaged by exposure to liquid. Do not allow liquids to come into contact with these windows.

## FUSES

This instrument does not contain fuses. If you suspect there is an electrical problem, please contact us for assistance.

# Corrosives, solvents, and pressurized gases

Many standard spectroscopy methods are based on the use of solvents. Others involve corrosive samples or pressurized samples in a gaseous state. All these sample materials can be measured using your instrument, but special precautions must be taken.

## VOLATILE SOLVENTS

If you use volatile solvents regularly, purging the instrument with dried, clean air or nitrogen is strongly recommended to create positive pressure inside the instrument. To avoid fire, explosion, and asphyxiation hazards, follow these guidelines when working with solvents.

- 
- Do not leave exposed solvent in the sample compartment for longer than necessary.
  - Work with the sample compartment cover open. For most infrared and near-infrared applications, it is not necessary to close the sample compartment cover.
  - Do not leave the solvents near the instrument.
  - If possible, purge the sample compartment.
  - Be sure that your work space is properly ventilated.

## DANGER

**Avoid fire and explosion.** The infrared source inside the instrument is an ignition source. If you are using volatile solvents, provide a fume hood or other active venting system that is free of spark and other ignition sources and prevents flammable vapors from collecting in the atmosphere surrounding the instrument.

These measures will help prolong the life of your instrument and will eliminate the possibility of spectral interference caused by volatile solvent vapors.

## WARNING

**Avoid hazard.** Materials such as hydrochloric acid, hydrofluoric acid and phosgene are highly toxic. If you regularly use solvents containing halogenated hydrocarbons, be sure your work area is properly ventilated.

## CORROSIVE SOLVENTS

Using corrosive solvents in the sample compartment may severely damage the system.

If you use corrosive solvents, such as those listed below, purging the instrument with dry, clean air or nitrogen is strongly recommended. Equipment damage due to failure to purge is not covered under the warranty. (If you have questions about this, please contact us.) The following is a list of commonly used halogenated solvents:

- Freon
- Methylene chloride
- Trichloroethylene
- Chloroform
- Carbon tetrachloride

---

**CONTACT US:**

[www.thermofisher.com](http://www.thermofisher.com)

Thermo Electron Scientific Instruments LLC

5225 Verona Rd.

Madison WI 53711 USA

PN 269-3648 00\_A

All rights reserved. Thermo Fisher Scientific Inc.

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