



ThermoFisher
S C I E N T I F I C

The world leader in serving science

DXR Data Collection Introduction



DXR User Training



Scope

- Best practices
 - Laser power
 - Aperture
- Automatic operation
 - Auto focus
 - Auto exposure
- Optimizing performance
 - Exposure time vs. number of exposures
 - Sampling
 - Microscope objectives
 - Smart Raman accessories



Best Practices

■ Laser power

- Use the maximum laser power whenever possible
- Situations requiring less laser power
 - Sample changes
 - Burning – signal disappears
 - Heating – signal or spectrum changes
 - Fluorescence

■ Aperture

- Use a slit aperture whenever possible
 - Pinhole and slit yield the same spectral resolution
 - Pinhole is required only for confocal operation
- Use a larger aperture if possible
 - Some loss of spectral resolution but it's not as bad as you may think
 - Try both aperture sizes and see if you lose any significant spectral resolution



Apertures


- Pinhole versus Slit

- Definition – aperture: opening through which light will pass
 - Pinhole – small, sharp edged opening, good depth of field
 - Slit – rectangular opening, large length versus narrow width
- Resolution – smaller opening, higher resolution, but loss of signal
- Confocal – use a pinhole to reject light from regions that are not in focal plane
- Applications
 - pinhole for microscopy work or small samples, for high resolution
 - slit for regular samples, bulk powders, liquids, for more signal

Apertures


■ Apertures

- 4 options
 - 2 pinholes and 2 slits
 - 25 and 50 micron
- Software selected
- Pinhole for confocal
- Slit for bulk




Effect of aperture size on spectral resolution

- 25um aperture spectrum superimposed on a 50um aperture spectrum (acetaminophen tablet)



DXR Spectral resolution comparison

- L-Cystine, 532 nm laser excitation



Auto exposure

Experiment Setup - C:\My Documents\omnic\VRParam\dxr_dc.exp

Collect | Bench | Quality | Advanced | Alignment | Mapping | Series

Estimated time for this collection: 00:02:01

Collect exposure time (sec): 1.0

Preview exposure time (sec): 1.0

Sample exposures: 2

Background exposures: 32

Final format: Shifted spectrum (cm-1)

Correction: Fluorescence

Cosmic ray threshold: Medium

Photobleach time (min): 0.0

Preview data collection

Auto exposure Desired S/N: 100

Maximum collect time (min): 2

File Handling

Save automatically Base name: []

Background Handling

Collect background before each sample

Maximum age for background: 1000 minutes


Use smart background

Experiment title:
DXR Data Collection

Experiment description:
Experiment file for DXR Data Collection training module

Help Open Save Save As OK Cancel

Autofocus





Optimizing performance

- Exposure time vs. number of exposures
 - Increasing the exposure time reduces the noise level
 - The signal level stays constant (counts per second)
 - Increasing the number of exposures reduces noise level
 - Random noise is reduced by signal averaging
- Both reduce noise; which is best?

Optimizing performance


- Exposure time vs. number of exposures
 - Increasing the exposure time reduces the noise level
 - The signal level stays constant (counts per second)
 - Increasing the number of exposures reduces noise level
 - Random noise is reduced by signal averaging
- Both reduce noise; which is best?

- Increasing exposure time reduces the noise level quicker than increasing the number of exposures

- For best results:
 - Leave # exposures = 2 (for cosmic ray rejection)
 - Increase the exposure time until CCD overflow


Comparison of exposure time and number of exposures

- Sugar, 532, 10x6-sec (red) vs. 2x30-sec exposure (blue)




Noise reduction without fluorescence

- Silicon, DXR Smart Raman, 532 nm laser excitation
- 60x1-sec (top, green)
- 30x2-sec (middle, blue)
- 2x30-sec (bottom, red)



Backgrounds

- Raman is an emission experiment so the concept of a background is much different than FT-IR which is (usually) an absorbance experiment
- A Raman background is a measure of the detector dark signal
 - Dark signal is the product of CCD temperature and exposure time
 - For the DXR, the detector temperature is constant
 - Measured with no light entering spectrograph (dark)
- Raman backgrounds are saved and reused until they expire
 - Background Handling control is on Collect tab of Experiment Setup



Background Handling

Collect background before each sample

Maximum age for background: minutes

Use smart background



Smart Backgrounds

- Patented technology that measures backgrounds at various exposure times from 1 to 60 seconds. A mathematical model is constructed which allows us to predict the background response for any exposure time
- Smart backgrounds
 - Eliminates need to collect backgrounds each time we change exposure time
 - Allows Auto Exposure capability based on *any* exposure time
 - Takes approximately 13 hours (overnight) to collect data for the model

Background handling controls

Experiment Setup - C:\My Documents\omnic\VRParam\dxr_dc.exp

Collect | Bench | Quality | Advanced | Alignment | Mapping | Series

Estimated time for this collection: 00:02:01

Collect exposure time (sec): 1.0

Preview exposure time (sec): 1.0

Sample exposures: 2

Background exposures: 32

Final format: Shifted spectrum (cm-1)

Correction: Fluorescence

Cosmic ray threshold: Medium

Photobleach time (min): 0.0

Preview data collection

Auto exposure Desired S/N: 100

Maximum collect time (min): 2

File Handling

Save automatically Base name: []

Background Handling

Collect background before each sample

Maximum age for background: 1000 minutes

Use smart background

Experiment title: DXR Data Collection

Experiment description: Experiment file for DXR Data Collection training module

Help Open Save Save As OK Cancel


Other topics

- We have covered the primary features contributing to data collection quality
 - Exposure time
 - Number of exposures
 - Aperture
 - Laser power
 - Autofocus

- The rest of this module will explore a variety of topics related more to operation and convenience
 - Final format
 - Cosmic ray rejection
 - Photobleaching
 - Corrections
 - Preview collect


Final format

- Selects the format for the collected spectrum
- Shifted spectrum (cm^{-1})
 - Raman shift
 - Typical format for Raman spectra
 - *Use this setting*
- Raman spectrum
 - Unshifted X-axis in either nm or cm^{-1}
 - Raman spectra collected with different lasers will not overlap
- Photoluminescence (nm)
 - X-axis in nanometers, Y-axis in emission units




Corrections

- Selects a spectral operation to apply automatically after collect and before displaying result
- None
- Raman efficiency
 - Generally not needed – advanced technique
 - Applies v^4 factor to Y-axis intensity
- Fluorescence
 - Removes fluorescence background using a polynomial fit or a reference spectrum
 - *Recommended setting (order=5)*



Fluorescence Correction Example



Uncorrected spectrum


Spectrum after automated fluorescence correction applied

Cosmic ray threshold

- Cosmic rays are natural events detected by sensitive integrating detectors such as a CCD.
 - Random, unpredictable events – not a specific wavelength
 - Typically show up as spikes in Raman spectra or backgrounds
- Rejection
 - Requires two exposures
 - Exposures are compared to detect and repair events
 - Turn-off (Threshold = None) to collect just one exposure
- Sensitivity threshold
 - Low is most sensitive (detects low energy events)
 - *Medium is recommended setting*

Cosmic Ray


- Here is an example of a spectrum that has a cosmic ray



Calcite spectrum



Preview collect

- Preview is a live display that displays a Raman spectrum so you can adjust sample focus or placement
 - Bench tab
 - Preview collect window before Collect
- Preview exposure time (sec):
 - Applies only to preview
 - 0.1 – 10 sec
- Preview data collection
 - Select (check) this option to automatically display a Preview Collect window before every data collection



Photobleaching

- Photobleaching exposes the sample to laser for a specified amount of time before starting data collection
- Reduces fluorescence by “bleaching” impurities causing the fluorescence





Alignment and calibration

- Two important but different processes for optimal performance
- Alignment is the process of bringing the laser beam and spectrograph sampling point into agreement
 - For the DXR Raman microscope, a third constraint is the visual crosshair which designates the sampling point
 - Alignment needed for best spatial resolution of microscope
- Calibration assures the accuracy of the spectrum wavelength (x-) axis and intensity (y-) axis
 - Software operations using known reference materials
 - Neon and polystyrene for the x-axis
 - White light for the y-axis






Calibrate Instrument

- The system is aligned. Next step: Calibration
- Calibration assures the accuracy of the spectrum wavelength (x-) axis and intensity (y-) axis
 - Software operations using known reference materials
 - Neon and polystyrene for the x-axis
 - White light for the y-axis
- Choose Collect > Calibrate Instrument from the OMNIC menu
- Place the alignment tool on the microscope stage and focus on the pinhole -or-
- Place the calibration accessory into the Smart Raman

Collect Sample

- Choose Collect > Collect Sample from the OMNIC menu



Add collect result to window

- The Confirmation dialog is displayed when the collect is finished
- The spectrum is still in the [Collect Sample] window
- If the result looks OK, press the Enter key or click the Yes button
- The spectrum is added to the active spectrum window

