

Synthetic Fibers by Microscope FT-IR Spectral Library

Index of Compounds

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
ACETATES					
2	AVTEX AVTEX	A0073	17-21 MICRONS	MULTILOBED, BRIGHT, WHITE	5990
3	AVTEX AVTEX	A0074	16-25 MICRONS	MULTILOBED, BRIGHT, WHITE	6090
4	AVTEX AVTEX	A0075	29-35 MICRONS	MULTILOBED, BRIGHT, WHITE	850
5	AVTEX AVTEX	A0076	14-22 MICRONS	MULTILOBED, DULL, WHITE	6941
6	AVTEX AVTEX	A0077	19-23 MICRONS	MULTILOBED, DULL, WHITE	6094
7	AVTEX AVTEX	A0078	17-26 MICRONS	MULTILOBED, DULL, WHITE	6550
8	AVTEX AVTEX	F0079	30-40 MICRONS	MULTILOBED, ----, BROWN	85-006-6
18	CELANESE ARNEL	A0300	21.6 MICRONS	IRREGULAR, BRIGHT, WHITE	SELVEDGE-3
19	CELANESE CELANESE	A0418	----	DOGBONE, BRIGHT IN WHITE LIGHT, WHITE	T-71
9	CELANESE TYPE REGULAR	A0287	21.3 MICRONS	IRREGULAR, BRIGHT, WHITE	REGULAR-1
10	CELANESE TYPE REGULAR	A0288	26.1 MICRONS	IRREGULAR, BRIGHT, WHITE	REGULAR-3
14	CELANESE TYPE REGULAR	A0295	19.0 MICRONS	IRREGULAR, BRIGHT, WHITE	REGULAR-5
16	CELANESE TYPE SELVEDGE	A0298	21.3 MICRONS	IRREGULAR, BRIGHT, WHITE	SELVEDGE-1
17	CELANESE TYPE SELVEDGE	A0299	21.3 MICRONS	IRREGULAR, DULL, WHITE	SELVEDGE-2
11	CELANESE TYPE D-151	A0289	21.3 MICRONS	IRREGULAR, BRIGHT, WHITE	D-151
12	CELANESE TYPE D-155	A0290	21.3 MICRONS	IRREGULAR, BRIGHT, WHITE	D-155
13	CELANESE TYPE D-162	A0292	21.3 MICRONS	IRREGULAR, BRIGHT, WHITE	D-162
15	CELANESE TYPE TOW	A0297	31.1 MICRONS	'Y' SHAPED TRI., DULL, WHITE	TOW-3
1	EASTMAN CHEMICAL ESTRON	A0067	26 MICRONS	MULTILOBED, BRIGHT, WHITE	REGULAR
348	WHITE ESTRON	A0066	26 MICRONS	MULTILOBED, DULL, WHITE	REGULAR
35	AMERICAN CYANAMID CRESLAN	A0160	24.5 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 61
ACRYLICS					
340	AMERICAN CYANAMID CRESLAN	A0159	24.5 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 68
21	BADISCHE (BASF) NO TRADE NAME	A0105	17.4 MICRONS	ROUND, SEMI-DULL, WHITE	A201
22	BADISCHE (BASF) NO TRADE NAME	A0106	16-19 MICRONS	ROUND, BRIGHT, WHITE	A201
23	BADISCHE (BASF) NO TRADE NAME	C0107	18-23 MICRONS	ROUND, SEMI-BRIGHT BLACK, BLACK	A405
24	BADISCHE (BASF) NO TRADE NAME	F0111	19.0 MICRONS	ROUND, BRIGHT, BROWN	A302
25	BADISCHE (BASF) NO TRADE NAME	A0110	15-24 MICRONS	ROUND, SEMI-BRIGHT, WHITE	A507
354	BADISCHE (BASF) NO TRADE NAME	L0109	19.0 MICRONS	ROUND, SEMI-DULL, BLUE	A451
39	E. I. DUPONT NOMEX	A0030	10-20 MICRONS	PEANUT, BRIGHT, WHITE	TYPE 450
349	E. I. DUPONT ORLON	A0011	18-33 MICRONS	MUSHROOM, SEMI-DULL, WHITE	TYPE 21
350	E. I. DUPONT ORLON	A0012	30-34 MICRONS (WDST)	DOGBONE, SEMI-DULL, WHITE	TYPE 42

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
351	E. I. DUPONT ORLON	A0028	18-26 MICRONS	MUSHROOM, SEMI-DULL, WHITE	TYPE 78
352	E. I. DUPONT ORLON	A0031	16-24 MICRONS	TRILOBAL, BRIGHT, WHITE	TYPE TR01
20	MONSANTO ACRILAN	A0100	24.6 MICRONS	---•---, SEMI-DULL, WHITE	TYPE 16
26	MONSANTO ACRILAN	A0113	24.6 MICRONS	BEAN, BRIGHT, WHITE	TYPE 16
27	MONSANTO ACRILAN	A0115	16-21 MICRONS	BEAN, BRIGHT, WHITE	TYPE 16
28	MONSANTO ACRILAN	A0114	19 MICRONS	BEAN, BRIGHT, WHITE	TYPE 16
29	MONSANTO ACRILAN	A0115	16-21 MICRONS	BEAN, BRIGHT, WHITE	TYPE 16
30	MONSANTO ACRILAN	A0117	12-17 MICRONS	BEAN, BRIGHT, WHITE	TYPE 16
31	MONSANTO ACRILAN	C0120	18-28 MICRONS	BEAN, BRIGHT BLACK, BLACK	TYPE 37
32	MONSANTO ACRILAN	L0121	18-28 MICRONS	BEAN, BRIGHT, BLUE	TYPE 39
33	MONSANTO ACRILAN	C0122	37-45 MICRONS	BEAN, BRIGHT BLACK, BLACK	TYPE 99
34	MONSANTO ACRILAN	A0124	42.6 MICRONS	BEAN, BRIGHT, WHITE	TYPE 90
337	MONSANTO ACRILAN	A0118	16-26 MICRONS	ROUND, BRIGHT, WHITE	TYPE X-24
338	MONSANTO ACRILAN	A0119	10-24 MICRONS	BEAN, BRIGHT, WHITE	TYPE 57
339	MONSANTO ACRILAN	A0123	34-50 MICRONS	BEAN-OVAL, BRIGHT, WHITE	TYPE 94
353	MONSANTO ACRILAN	A0101	---•---	---•---, ---•---, WHITE	TYPE 16
ARAMIDS					
36	E. I. DUPONT KEVLAR 29	H0009	11-13 MICRONS	ROUND, BRIGHT, YELLOW	TYPE 964
38	E. I. DUPONT KEVLAR 49	H0029	11-13 MICRONS	ROUND, BRIGHT, YELLOW	TYPE 966
37	E. I. DUPONT NOMEX	A0014	10-20 MICRONS	PEANUT, BRIGHT, WHITE	200-100-TYPE 431
355	E. I. DUPONT NOMEX	A0030	---•---	---•---, ---•---, WHITE	TYPE 450
FLOUROCARBONS					
346	ALBANY INTERNATIONAL K4-9	A0465	---•---	ROUND, ---•---, WHITE	005 K4-9
347	ALBANY INTERNATIONAL T2-3	D0466	---•---	ROUND, ---•---, PINK	008 T2-3
334	E. I. DUPONT TEFLON	A0023	27-30 MICRONS	ROUND, BRIGHT, WHITE	1350-180-0
335	E. I. DUPONT TEFLON	F0024	27-30 MICRONS	ROUND, DULL, BROWN	200-30-0
345	W. L. GORE & ASSOC. GORE-TEX	A0431	---•---	---•---, ---•---, WHITE	---•---
MODACRYLICS					
40	EASTMAN CHEMICAL VEREL	A0068	14-56 MICRONS	DOGBONE, DULL, WHITE	TYPE 163
41	MONSANTO ACRILAN	A0125	16-20 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE SEF
42	MONSANTO SEF	A0381	---•---	ROUND, SEMI-DULL, WHITE	S-06
43	MONSANTO SEF	A0382	---•---	ROUND, SEMI-DULL, WHITE	S-06
NYLONS					
372	ALBANY INTERNATIONAL N14-8	A0463	---•---	ROUND, ---•---, WHITE	008 N14-8

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
123	ALBANY INTERNATIONAL N4-9	A0464	----	ROUND, ----, WHITE	008 N4-9
87	ALLIED ANSO	A0218	45.2 MICRONS	TRI-SYMMETRIC, BRIGHT, WHITE	1183-70 7B39
74	ALLIED ANSO IV	A0198	45.2 MICRONS	TRI-SYMMETRIC, BRIGHT, WHITE	1185-70-00 7B66
77	ALLIED ANSO IV	A0202	45.1 MICRONS	TRI-ASYMMETRIC, SEMI-DULL, WHITE	1180-70-00 6N67
78	ALLIED ANSO IV	A0203	46.0 MICRONS	TRI-ASYMMETRIC, MID-DULL, WHITE	1225-70-00 6X67
79	ALLIED ANSO IV	A0204	45.3 MICRONS	TRI-ASYMMETRIC, BRIGHT, WHITE	1784-105-00 6B68
86	ALLIED ANSO IV	A0217	45.3 MICRONS	SYMMETRIC 'Y', BRIGHT, WHITE	NO CODE
75	ALLIED ANSO IV AM HALOFRES	A0199	57.0 MICRONS	TRI-SYMMETRIC, MID-DULL PERFECT LUSTER, WHITE	2582-96-00 7J31
73	ALLIED ANSO IV AM	A0197	45.2 MICRONS	TRI-SYMMETRIC, BRIGHT, WHITE	1183-70-00 7B68
76	ALLIED ANSO IV AM	A0200	56.8 MICRONS	TRI-SYMMETRIC, MID-DULL PERFECT LUSTER, WHITE	3840-144-00 7K51
83	ALLIED CAPROLAN	H0211	27.3 MICRONS	ROUND, BRIGHT, YELLOW	1260-204.3Z N62-24
84	ALLIED CAPROLAN	H0213	----	----, ----, YELLOW	1260-204.5Z N62-24
360	ALLIED CAPROLAN	H0205	27.3 MICRONS	ROUND, BRIGHT, YELLOW	1260-204.5Z N14-20
362	ALLIED CAPROLAN	H0206	----	ROUND, BRIGHT, YELLOW	N50-20
363	ALLIED CAPROLAN	H0212	27.3 MICRONS	ROUND, BRIGHT, YELLOW	1260-204.3Z N62-24
364	ALLIED CAPROLAN	A0269	27.3 MICRONS	ROUND, BRIGHT, WHITE	N67-24 BHS
365	ALLIED CAPROLAN	A0271	13.5 MICRONS	ROUND, SEMI-DULL, WHITE	M-1460
366	ALLIED CAPROLAN	A0272	22.0 MICRONS	DELTA, BRIGHT, WHITE	M-1470
367	ALLIED CAPROLAN	A0273	26.6 MICRONS	ROUND, MID-DULL, WHITE	M-791
368	ALLIED CAPROLAN	A0275	20.1 MICRONS	DELTA, BRIGHT, WHITE	M-425
369	ALLIED CAPROLAN	A0277	15.8 MICRONS	ROUND, MID-DULL, WHITE	M-737
85	ALLIED LUROL	A0215	45.3 MICRONS	SYMMETRIC 'Y', BRIGHT, WHITE	28326
80	ALLIED NO TRADE NAME	A0207	27.3 MICRONS	ROUND, BRIGHT, WHITE	840-136.3Z N50-20
81	ALLIED NO TRADE NAME	A0209	31.6 MICRONS	ROUND, BRIGHT, WHITE	1680-204.5Z N50-20
82	ALLIED NO TRADE NAME	A0210	27.3 MICRONS	ROUND, BRIGHT, WHITE	1260-204.3Z N52-20
361	ALLIED NO TRADE NAME	A0208	27.3 MICRONS	ROUND, BRIGHT, WHITE	1260-204.3Z N50-20
61	BADISCHE (BASF) ----	A0162	20.1 MICRONS	----, ----, WHITE	1832
62	BADISCHE (BASF) ----	A0163	----	----, ----, WHITE	----
63	BADISCHE (BASF) ----	A0164	----	----, ----, WHITE	----
64	BADISCHE (BASF) NO TRADE NAME	A0165	60.4 MICRONS	PENTAGONAL, MID-DULL, WHITE	T255
89	BADISCHE (BASF) ZEFTRON	E0248	51.6 MICRONS	TRI-SYMMETRIC, SEMI-DULL, RED	TYPE W-132
65	BADISCHE (BASF) ZEFTRON 500	B0166	46.5 MICRONS	ROUND, SUBDUED, GRAY	T23A-403
66	BADISCHE (BASF) ZEFTRON 500 ZX	0167	45.9 MICRONS	PENTAGONAL, SEMI-DULL, MULTICOLOR	T31B-503
88	BADISCHE ZEFTRON 500 ZX	F0247	45.9 MICRONS	PENTALOBAL, SEMI-DULL, BROWN	TYPE H-110
112	BERKLEY AND CO ----	A0408	----	ROUND, CLEAR, WHITE	----

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
117	BERKLEY AND CO DEWFLEX	A0410	---•---	ROUND, CLEAR, WHITE	---•---
114	BERKLEY AND CO TRILENE	A0409	---•---	ROUND, OPAQUE, WHITE	---•---
119	CAMAC CORPORATION CAMALON	B0429	---•---	'Y', SEMI-DULL, GRAY	NONE
120	CAMAC CORPORATION CAMALON	F0430	---•---	ROUND, SEMI-DULL, BROWN	NONE
45	E. I. DUPONT ANTRON	A0008	12 MICRONS	TRILOBAL, BRIGHT, WHITE	140-34-TYPE 860
47	E. I. DUPONT ANTRON	A0018	44 MICRONS	TETRALOBAL, SLIGHTLY DELUSTERED, WHITE	1225-80-TYPE 857
57	E. I. DUPONT ANTRON III	A0020	47 MICRONS	TETRALOBAL, BRIGHT, WHITE	2470-131-TYPE 756A
48	E. I. DUPONT NO TRADE NAME	A0019	51-65 MICRONS	TRILOBAL, BRIGHT, WHITE	1350-68-TYPE 746
49	E. I. DUPONT NO TRADE NAME	A0103	17.4 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 420
50	E. I. DUPONT NO TRADE NAME	A0104	19 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 200
44	E. I. DUPONT QIANA	A0007	18 MICRONS	TRILOBAL, SLIGHTLY DELUSTERED, WHITE	60-36-TYPE 472
72	ENKA (BASF) ANKALON	A0187	51.6 MICRONS	TRIANGULAR, SEMI-DULL, WHITE	NO CODE
68	ENKA (BASF) CREPE SET	A0175	60.2 MICRONS	OVAL, DULL, WHITE	8577
69	ENKA (BASF) ENKALON	A0179	47.7 MICRONS	TRILOBAL, BRIGHT, WHITE	8914
70	ENKA (BASF) ENKALON	A0180	46.2 MICRONS	TRILOBAL, BRIGHT, WHITE	8380
71	ENKA (BASF) ENKALON	A0186	46.7 MICRONS	TRILOBAL, BRIGHT, WHITE	NO CODE
90	ENKA (BASF) ENKALON	A0261	44.1 MICRONS	---•---, SEMI-BRIGHT, WHITE	3817
91	ENKA (BASF) ENKALON AP	A0262	51.4 MICRONS	TRILOBAL, MIDLUSTER, WHITE	3435
92	ENKA (BASF) ENKALON AP	A0263	51.6 MICRONS	---•---, BRIGHT, WHITE	---•---
358	ENKA (BASF) ENKALURE	A0177	16.3 MICRONS	IRREGULAR TRI, BRIGHT, WHITE	8378
359	ENKA (BASF) ENKALURE	A0178	22.4 MICRONS	TRIANGULAR, BRIGHT, WHITE	8770
67	ENKA (BASF) ENKASHEER	A0172	42.6 MICRONS	DISTORTED ROUND, SEMI-DULL, WHITE	8463
357	ENKA (BASF) NO TRADE NAME	A0174	15.6 MICRONS	---•---, SEMI-DULL, WHITE	8857
356	ENKA (BASF) SOFTALON	A0173	22.0 MICRONS	DISTORTED ROUND, SEMI-DULL, WHITE	7614
118	KAYSER-ROTH NONE	A0422	---•---	ROUND, BRIGHT, WHITE	---•---
371	MONOFILAMENTS MONOFUSE A/S	A0307	228.6 MICRONS	ROUND, CLEAR, WHITE	#40
370	MONOFILAMENTS STEAM-AWAY	A0302	31.7 MICRONS	ROUND, CLEAR, WHITE	150D
51	MONSANTO ACTION WEAR	A0127	15.8 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE L80
55	MONSANTO NO TRADE NAME	A0131	46.7 MICRONS	TRILOBAL, BRIGHT, WHITE	T-1878-CBS
52	MONSANTO ULTRON	A0128	31-34 MICRONS	TRILOBAL, BRIGHT, WHITE	T-1091-CES
53	MONSANTO ULTRON	A0129	42.6 MICRONS	TRILOBAL, BRIGHT, WHITE	T-1551-CBS
59	MONSANTO ULTRON	A0134	42.3 MICRONS	TRILOBAL, BRIGHT, WHITE	T-1070-DBE
60	MONSANTO ULTRON	A0136	42.6 MICRONS	TRILOBAL, BRIGHT, WHITE	T-1275-CBE
98	MONSANTO ULTRON	A0383	---•---	TRILOBAL, MID-DULL, WHITE	T-1233
99	MONSANTO ULTRON	A0384	---•---	TRILOBAL, MID-DULL, WHITE	T-1633

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
100	MONSANTO ULTRON	A0385	---	TRILOBAL, MID-DULL, WHITE	T-1775
101	MONSANTO ULTRON	A0386	---	TRILOBAL, SUPER BRIGHT, WHITE	T-1891
102	MONSANTO ULTRON	A0387	---	TRILOBAL, MID-DULL, WHITE	T-1963
103	MONSANTO ULTRON	A0388	---	TRILOBAL, MID-DULL, WHITE	T-1983
104	MONSANTO ULTRON	A0389	---	TRILOBAL, MID-DULL, WHITE	T-1985
105	MONSANTO ULTRON	A0390	---	TRILOBAL, BRIGHT, WHITE	T-1988
106	MONSANTO ULTRON	A0391	---	TRILOBAL, SEMI-DULL, WHITE	T-1993
107	MONSANTO ULTRON	A0392	---	TRILOBAL, SEMI-DULL, WHITE	T-2089
108	MONSANTO ULTRON	A0393	---	TRILOBAL, SEMI-DULL, WHITE	T-2090
109	MONSANTO ULTRON	A0394	---	TRILOBAL, SEMI-DULL, WHITE	T-2093
110	MONSANTO ULTRON	A0395	---	TRILOBAL, SUPER BRIGHT, WHITE	1340/68 CEJ
54	MONSANTO ULTRON Z	A0130	42.3 MICRONS	TRILOBAL, SEMI-DULL, WHITE	T-1595-LZS
58	MONSANTO ULTRON Z 3D	A0133	49.2 MICRONS	TRILOBAL, SEMI-DULL, WHITE	T-2095-LZS
115	SHAKESPEARE CO. ---	H0401	---	ROUND, OPAQUE SOLUTION DYED, YELLOW	MN7
56	MONSANTO ULTRON 3D	A0132	49.2 MICRONS	TRILOBAL, SEMI-DULL, WHITE	T-2000-LES
113	SHAKESPEARE CO. ---	A0399	---	ROUND, CLEAR, WHITE	WN-18
116	SHAKESPEARE CO. ---	A0403	---	ROUND, ---, WHITE	FN-29
111	SHAKESPEARE CO. WONDER THREAD	A0402	---	ROUND, CLEAR, WHITE	SN38
95	STAR FIBERS INC. NONE	A0356	---	MULTILOBAL, BRIGHT, WHITE	2900A
96	STAR FIBERS INC. NONE	A0357	---	MULTILOBAL, BRIGHT, WHITE	1982-A
97	STAR FIBERS INC. NONE	A0358	---	MULTILOBAL, SEMI-DULL, WHITE	1983
94	TOLARAM FIBERSINC ---	A0330	---	ROUND, SEMI-DULL, WHITE	---
121	WELLMAN, INC. WELLON	A0442	---	TRILOBAL, SEMI-DULL, WHITE	TYPE 180T
122	WELLMAN, INC. WELLON	A0444	---	TRILOBAL, BRIGHT, WHITE	TYPE 180BT

OLEFINS

170	ALADDIN MILLS, INC. XYLON	0440	---	DELTA, SEMI-DULL, MULTICOLOR	---
137	AMOCO ---	H0258	---	FLAT, DULL, YELLOW	1273
138	AMOCO ---	H0259	---	FLAT, BRIGHT, YELLOW	1273
133	AMOCO AMOCO STAPLE	E0254	19.0 MICRONS	DELTA, DULL, RED	TYPE 013
136	AMOCO AMOCO OLEFIN	E0257	46.7 MICRONS	DELTA, DULL, RED	TYPE 010
134	AMOCO AMOCO CF YARN	F0255	17.7 MICRONS	ROUND, BRIGHT, BROWN	TYPE 030
135	AMOCO AMOCO CF YARN	H0256	22.5 MICRONS	DELTA, BRIGHT, YELLOW	TYPE 010
373	AMOCO SYNERA	L0253	14.8 MICRONS	ROUND, BRIGHT, BLUE	TYPE 010
144	AVTEX ---	A0313	15.6 MICRONS	ROUND, BRIGHT, WHITE	PNE 7026
168	BEAULIEU OF AMERICA NONE	F0437	---	TRILOBAL, ---, BROWN	NONE

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
160	BERKLEY AND CO BERKLEY POLYPROPYLENE	A0406	----	ROUND, CLEAR, WHITE	----
161	BERKLEY AND CO BERKLEY POLYETHYLENE	A0407	----	ROUND, OPAQUE, WHITE	----
157	CHASE BAG COMPANY POLYSAX-X	E0398	----	----, OPAQUE IN NATURAL, RED	----
124	CHEVRON CHEM. (VECTRA) POLYLOOM	K0032	VARIABLE	FLAT, BRIGHT, GREEN	5700 DENIER
125	CHEVRON CHEM. (VECTRA) POLYLOOM	K0033	VARIABLE	FLAT, BRIGHT, GREEN	3800 DENIER
126	CHEVRON CHEM. (VECTRA) VECTRA	A0034	25-38 MICRONS	ROUND, SEMI-DULL, WHITE	400 DENIER
127	CHEVRON CHEM. (VECTRA) VECTRA	A0035	30-40 MICRONS	TRILOBAL, SLIGHTLY DELUSTERED, WHITE	300 DENIER
128	CHEVRON CHEM. (VECTRA) VECTRA	A0036	30-38 MICRONS	ROUND-MODIFIED, SLIGHTLY DELUSTERED, WHITE	7 DENIER
158	CONCORDE FIBERS, INC. NONE	H0404	----	DELTA, BRIGHT, YELLOW	2600/108
159	CONCORDE FIBERS, INC. NONE	L0405	----	ROUND, BRIGHT, BLUE	2400/160
164	CUMBERLAND MILLS, INC. NONE	B0415	----	DELTA, DULL, GRAY	2400
146	ELIZABETH WEBBING NONE	E0370	----	----, DULL IN UNDYED STATE, RED	FORTILENE 3807
172	ENSIGN- BICKFORD ----	A0448	----	----, OPAQUE, WHITE	----
173	ENSIGN- BICKFORD SUPERBULK CABLE FILL	L0449	----	----, OPAQUE, BLUE	----
188	EXXON CHEMICAL ----	H0478	----	RECTANGULAR, ----, YELLOW	CBW WARP
189	EXXON CHEMICAL ----	C0479	----	RECTANGULAR, ----, BLACK	----
190	EXXON CHEMICAL ----	A0480	----	RECTANGULAR, ----, WHITE	IBC 2000
169	FIBERS SOUTH, INC. ----	A0438	----	----, ----, WHITE	----
155	FIBRON CORP. FIBRILON	A0396	----	----, ----, WHITE	NONE
93	HERCULES ----	A0264	----	----, ----, WHITE	T-116
139	HERCULES HERCULON	A0265	19.0 MICRONS	ROUND, BRIGHT, WHITE	T-101
140	HERCULES HERCULON	A0266	22.4 MICRONS	DELTA, BRIGHT, WHITE	T-335
141	HERCULES HERCULON	A0267	26.9 MICRONS	ROUND, BRIGHT, WHITE	T-303
142	HERCULES HERCULON	A0268	46 MICRONS	DELTA, BRIGHT, WHITE	T-701
166	INTEGRATED PRODUCTS ----	A0427	----	ROUND, ----, WHITE	----
156	J & J INDUSTRIES ----	N0397	----	----, ----, PURPLE	NONE
187	JAMES RIVER CORP. ----	A0471	----	ROUND, ----, WHITE	L2
175	LORETEX CORP. LORETEX	A0451	----	----, ----, WHITE	NONE
176	LORETEX CORP. LORETEX	A0454	----	----, ----, WHITE	4 X 5 WHITE
177	LORETEX CORP. LORETEX	A0455	----	----, ----, WHITE	NONE
178	LORETEX CORP. LORETEX	A0456	----	----, ----, WHITE	4 X D2 WHITE

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
179	LORETEX CORP. LORETEX	A0457	----	-----, -----, WHITE	8 X 8 WHITE
180	LORETEX CORP. LORETEX	L0458	----	-----, -----, BLUE	10 X 10 BLUE
181	LORETEX CORP. LORETEX	H0459	----	-----, -----, YELLOW	NONE
182	LORETEX CORP. LORETEX	A0461	----	-----, -----, WHITE	4 X D2
183	LORETEX CORP. LORETEX	H0460	----	-----, -----, YELLOW	12 X D6
375	LORETEX CORP. LORETEX	C0452	----	-----, -----, BLACK	NONE
376	LORETEX CORP. LORETEX	C0453	----	-----, -----, BLACK	12 X 12 BLACK
162	MIRAFI INC. MIRAFI	A0412	----	-----, -----, WHITE	MIRADRAIN 6000
374	MIRAFI INC. MIRAFI	C0413	----	-----, -----, BLACK	700X
143	MONOFILAMENTS -----	A0305	192 MICRONS	ROUND, CLEAR, WHITE	375D
171	OSTERNECK CO. -----	B0441	----	-----, -----, GRAY	NONE
129	PHILLIPS MARVSS	F0249	18.4 MICRONS	DELTA, HIGH LUSTER, BROWN	G14K09
130	PHILLIPS MARVSS	F0250	29.1 MICRONS	MODIFIED DELTA, REGULAR, BROWN	G14F69
131	PHILLIPS MARVSS	H0251	46.7 MICRONS	ROUND, REGULAR, YELLOW	G14J89
132	PHILLIPS MARVSS	H0252	85.2 MICRONS	ROUND, REGULAR, YELLOW	G14J63
163	POLYLOOM CORP. OF AMERICA POLYLOOM	K0414	----	-----, -----, GREEN	-----
167	SHUFORD MILLS SHUR-TIE	A0436	----	-----, -----, WHITE	S-2100
186	SOLUTION FIBERS, INC. LUXENE	F0470	----	DELTA, -----, BROWN	NONE
184	SUNBURST YARNS, INC. NONE	F0468	----	DELTA, -----, BROWN	1326
185	SUNBURST YARNS, INC. NONE	E0469	----	ROUND, -----, RED	863AE
174	SYNTHETIC INDUSTRIES -----	A0450	----	ROUND, DULL, WHITE	STYLE 407
145	WALTRICH PLASTIC CP NONE	A0344	----	-----, -----, WHITE	-----
165	WAYN-TEX INC. -----	F0419	----	FLAT, ROUND EDGES, -----, BROWN	-----
147	WELLINGTON PURITAN -----	F0372	----	-----, -----, BROWN	109 DIE-DELTA
148	WELLINGTON PURITAN -----	F0374	----	-----, -----, BROWN	DERBY ROPE
149	WELLINGTON PURITAN -----	C0375	----	-----, -----, BLACK	-----
150	WELLINGTON PURITAN -----	H0376	----	-----, -----, YELLOW	M1501
152	WELLINGTON PURITAN -----	C0378	----	-----, -----, BLACK	-----
154	WELLINGTON PURITAN -----	G0380	----	-----, -----, ORANGE	-----
153	WELLINGTON PURITAN CO-POLY	A0379	----	-----, -----, WHITE	-----
151	WELLINGTON PURITAN UNMANILA	F0377	----	-----, -----, BROWN	-----

POLYESTERS

239	ALLIED PET	A0188	29.4 MICRONS	ROUND, BRIGHT, WHITE	500-70-00-70 N25-5
-----	---------------	-------	--------------	----------------------	--------------------

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
240	ALLIED PET	A0189	38.1 MICRONS	ROUND, BRIGHT, WHITE	840-70-00-70/N43-30
241	ALLIED PET	L0190	25.1 MICRONS	ROUND, BRIGHT, BLUE	1000-192-00-70 N43-30
242	ALLIED PET	A0191	25.1 MICRONS	ROUND, BRIGHT, WHITE	1000-192-00-70 N44-100
243	ALLIED PET	A0192	25.1 MICRONS	ROUND, BRIGHT, WHITE	1000-192-00-72 N43-30
244	ALLIED PET	K0193	25.1 MICRONS	ROUND, BRIGHT, GREEN	1000-192-00-73 N51-33
245	ALLIED PET	A0195	28.6 MICRONS	ROUND, BRIGHT, WHITE	1300-192-00-73 N32-32
246	ALLIED PET	A0196	29.0 MICRONS	ROUND, BRIGHT, WHITE	1600-230-00-70 N23-30
228	AVTEX AVLIN	A0083	12 MICRONS	ROUND, SEMI-DULL, WHITE	7605
229	AVTEX AVLIN	A0084	18-24 MICRONS	TRILOBAL, SEMI-DULL, WHITE	9613
230	AVTEX AVLIN	A0086	17-18 MICRONS	ROUND, SEMI-DULL, WHITE	9308
231	AVTEX AVLIN	A0087	16-26 MICRONS	TRILOBAL, SEMI-DULL, WHITE	8666
272	AVTEX AVLIN	A0314	19.0 MICRONS	TRILOBAL, BRIGHT, WHITE	8438
273	AVTEX AVLIN	A0315	19.0 MICRONS	ROUND, BRIGHT, WHITE	9308
274	AVTEX AVLIN	A0316	---•---	TRILOBAL, BRIGHT, WHITE	7814
262	CELANESE FORTREL	A0278	12.0 MICRONS	ROUND, SEMI-DULL, WHITE	T-367
263	CELANESE FORTREL	A0279	13.5 MICRONS	ROUND, SEMI-DULL, WHITE	T-300
264	CELANESE FORTREL	A0280	13.5 MICRONS	ROUND, SEMI-DULL, WHITE	T-310
265	CELANESE FORTREL	A0281	12.7 MICRONS	TRILOBAL, BRIGHT, WHITE	T-702
266	CELANESE FORTREL	A0282	18.0 MICRONS	TRILOBAL, SEMI-DULL, WHITE	T-708
267	CELANESE FORTREL	A0283	15.6 MICRONS	ROUND, SEMI-DULL, WHITE	T-728-1
268	CELANESE FORTREL	A0284	15.6 MICRONS	ROUND, DULL, WHITE	T-728-2
269	CELANESE FORTREL	A0285	23.4 MICRONS	ROUND, BRIGHT, WHITE	T-700
294	CELANESE FORTREL	A0416	---•---	TRILOBAL, ---•---, WHITE	T-427
295	CELANESE FORTREL	A0417	---•---	ROUND, EXTRA WHITE, WHITE	T-472 OR T-478
291	COLOR-FI INC. NONE	A0355	---•---	MODIFIED TRI, BRIGHT, WHITE	15-101-7
292	COLOR-FI INC. NONE	A0369	---•---	ROUND, SEMI-DULL, WHITE	15-108-2R
336	E. I. DUPONT CANTRECE	A0025	30 MICRONS	ROUND, SLIGHTLY DELUSTERED, WHITE	45-7-TYPE 882
191	E. I. DUPONT DACRON	A0001	13-16 MICRONS	ROUND, DULL, WHITE	70-34-TYPE 57
192	E. I. DUPONT DACRON	A0003	24-40 MICRONS	DOUGHNUT, SLIGHTLY DELUSTERED, WHITE	HOLLOWFIL-TYPE 808
193	E. I. DUPONT DACRON	A0005	15 MICRONS	ROUND, BRIGHT, WHITE	70-34-TYPE 55
194	E. I. DUPONT DACRON	C0006	21 MICRONS	ROUND, BLACK-CARBON PIGMENT, BLACK	150-34-TYPE 55A
195	E. I. DUPONT DACRON	A0017	20-22 MICRONS	ROUND, SEMI-DULL, WHITE	150-34-TYPE 56
196	E. I. DUPONT DACRON	A0010	15-19 MICRONS	TRILOBAL, SEMI-DULL, WHITE	70-34-TYPE 62
197	E. I. DUPONT DACRON	A0021	18-23 MICRONS	ROUND, SLIGHTLY DELUSTERED, WHITE	150-34-TYPE 232
199	E. I. DUPONT DACRON	A0026	24-30 MICRONS	DOUGHNUT, SLIGHTLY DELUSTERED, WHITE	HOLLOWFIL II
200	E. I. DUPONT DACRON	A0027	22-30 MICRONS	DOUGHNUT, SLIGHTLY DELUSTERED, WHITE	HOLLOWFIL IIF

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
232	E. I. DUPONT DACRON	A0099	16.5 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 54
46	E. I. DUPONT NO TRADE NAME	A0016	17 MICRONS	ROUND, SEMI-DULL, WHITE	15-7-TYPE 288
220	EASTMAN CHEMICAL KODEL	A0061	27 MICRONS	ROUND, BRIGHT, WHITE	TYPE 411
221	EASTMAN CHEMICAL KODEL	A0063	23-30 MICRONS	ROUND TO OVAL, SEMI-DULL, WHITE	TYPE 435
222	EASTMAN CHEMICAL KODEL	A0064	43 MICRONS	TRILOBAL, BRIGHT, WHITE	TYPE 417
223	EASTMAN CHEMICAL KODEL	A0065	24 MICRONS	SWOLLEN RIBBON, BRIGHT, WHITE	TYPE 419
224	EASTMAN CHEMICAL KODEL	A0069	21-30 MICRONS	POLYGONAL, SEMI-DULL, WHITE	TYPE 30F
225	EASTMAN CHEMICAL KODEL	A0070	14-36 MICRONS	POLYGONAL, SEMI-DULL, WHITE	TYPE 32F
226	EASTMAN CHEMICAL KODEL	A0071	13-18 MICRONS	POLYGONAL, SEMI-DULL, WHITE	TYPE 40P
227	EASTMAN CHEMICAL KODEL	A0072	12-38 MICRONS	ROUNDED POLYGON, SEMI-DULL, WHITE	TYPE 30P
279	EASTMAN CHEMICAL KODEL	A0335	---•---	ROUND, BRIGHT, WHITE	TYPE 211
280	EASTMAN CHEMICAL KODEL	A0336	---•---	ROUND, SEMI-DULL, WHITE	TYPE 425
281	EASTMAN CHEMICAL KODEL	A0337	---•---	ROUND, SEMI-DULL, WHITE	TYPE 426
234	ENKA (BASF) ENKRON	A0181	16.3 MICRONS	---•---, BRIGHT, WHITE	7181
235	ENKA (BASF) ENKRON	A0182	17.4 MICRONS	DISTORTED ROUND, SEMI-DULL TEXTURED, WHITE	---•---
236	ENKA (BASF) ENKRON	A0183	19.4 MICRONS	TRIANGULAR, BRIGHT, WHITE	7151
237	ENKA (BASF) ENKRON	A0184	19.3 MICRONS	ROUND, SEMI-DULL, WHITE	7444
238	ENKA (BASF) ENKRON	A0185	24.6 MICRONS	DISTORTED ROUND, SEMI-DULL TEXTURED, WHITE	---•---
296	FIRESTONE FIBERS ---•---	A0420	---•---	ROUND, OPAQUE, WHITE	TYPE 005
297	FIRESTONE FIBERS ---•---	A0421	---•---	ROUND, OPAQUE, WHITE	TYPE 007
300	GOODYEAR TIRE & RUBBER 826P STD MOD TIRE YARN	A0439	23.1 MICRONS	CIRCULAR, OPAQUE, WHITE	826P
248	HOECHST HOECHST	A0233	16.5 MICRONS	ROUND, CLEAR, WHITE	210
250	HOECHST HOECHST	A0235	13.5 MICRONS	ROUND, SEMI-DULL, WHITE	290
251	HOECHST HOECHST	A0236	42.6 MICRONS	PENTALOBAL, CLEAR, WHITE	294
252	HOECHST HOECHST	A0237	42.6 MICRONS	PENTALOBAL, SEMI-DULL, WHITE	295
253	HOECHST HOECHST	A0238	42.6 MICRONS	TRILOBAL, CLEAR, WHITE	802
254	HOECHST HOECHST	A0239	42.6 MICRONS	TRILOBAL, CLEAR, WHITE	803
255	HOECHST HOECHST	A0240	46.7 MICRONS	PENTALOBAL, CLEAR, WHITE	814
256	HOECHST HOECHST	A0241	42.6 MICRONS	TRILOBAL, CLEAR, WHITE	817
257	HOECHST HOECHST	A0242	46.7 MICRONS	PENTALOBAL, SEMI-DULL, WHITE	885
286	HOECHST HOECHST	A0349	---•---	ROUND, SEMI-DULL, WHITE	TYPE 261
287	HOECHST HOECHST	A0350	---•---	ROUND, BRIGHT, WHITE	TYPE 290
288	HOECHST HOECHST	A0352	---•---	PENTALOBAL, SEMI-BRIGHT, WHITE	TYPE 826
289	HOECHST HOECHST	A0353	---•---	TRILOBAL, CLEAR, WHITE	TYPE 843
290	HOECHST HOECHST	A0354	---•---	PENTALOBAL, CLEAR, WHITE	TYPE 845

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
198	HOECHST TREVIRA	A0037	12 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 101
201	HOECHST TREVIRA	A0039	18 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 220
202	HOECHST TREVIRA	A0040	16 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 221
203	HOECHST TREVIRA	A0041	45 MICRONS	PENTALOBAL, SEMI-DULL, WHITE	TYPE 265
204	HOECHST TREVIRA	A0042	26 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 270
205	HOECHST TREVIRA	A0043	10 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 271
206	HOECHST TREVIRA	A0044	25 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 280
207	HOECHST TREVIRA	A0045	13 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 291
208	HOECHST TREVIRA	A0046	26 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 291
209	HOECHST TREVIRA	A0047	44 MICRONS	PENTALOBAL, BRIGHT, WHITE	TYPE 295
210	HOECHST TREVIRA	A0048	24 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 296
211	HOECHST TREVIRA	A0049	17 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE 350
212	HOECHST TREVIRA	A0051	24 MICRONS	PENTALOBAL, BRIGHT, WHITE	TYPE 815
213	HOECHST TREVIRA	A0052	40 MICRONS	TRILOBAL, BRIGHT, WHITE	TYPE 816
214	HOECHST TREVIRA	A0053	34 MICRONS	PENTALOBAL, BRIGHT, WHITE	TYPE 818
215	HOECHST TREVIRA	A0054	30-52 MCRONS	PENTALOBAL, BRIGHT, WHITE	TYPE 825
216	HOECHST TREVIRA	A0055	47 MICRONS	TRILOBAL, SEMI-DULL, WHITE	TYPE 873
217	HOECHST TREVIRA	A0058	28 MICRONS	POLYGONAL, SEMI-DULL, WHITE	TYPE 632
218	HOECHST TREVIRA	A0059	32 MICRONS	POLYGONAL, SEMI-DULL, WHITE	TYPE 632
219	HOECHST TREVIRA	A0060	22 MICRONS	MULT-POLYGONAL, SEMI-DULL, WHITE	TYPE 692
247	HOECHST TREVIRA	A0232	16.5 MICRONS	ROUND, SEMI-DULL, WHITE	141
249	HOECHST TREVIRA	A0234	19.0 MICRONS	ROUND, BRIGHT, WHITE	220
258	HOECHST TREVIRA	A0243	30.1 MICRONS	ROUND, SEMI-DULL, WHITE	618
259	HOECHST TREVIRA	A0244	23.8 MICRONS	OCTALOBAL, SEMI-DULL, WHITE	618
260	HOECHST TREVIRA	A0245	23.8 MICRONS	PENTALOBAL, SEMI-DULL, WHITE	625
261	HOECHST TREVIRA	A0246	23.8 MICRONS	PENTALOBAL, SEMI-DULL, WHITE	695
282	HOECHST TREVIRA	A0345	---•---	PENTALOBAL, SDOW, WHITE	TYPE 204
283	HOECHST TREVIRA	A0346	---•---	PENTALOBAL, SEMI-DULL, WHITE	TYPE 205
284	HOECHST TREVIRA	A0347	---•---	ROUND HOLLOW, SEMI-DULL, WHITE	TYPE 206
285	HOECHST TREVIRA	A0348	---•---	---•---, ---•---, WHITE	TYPE 207
275	METLON CORP. RETROGLO	B0329	---•---	---•---, ---•---, GRAY	NONE
270	MONOFILAMENTS NO TRADE NAME	A0303	76.2 MICRONS	ROUND, CLEAR, WHITE	#5
271	MONOFILAMENTS NO TRADE NAME	A0304	355 MICRONS	ROUND, CLEAR, WHITE	1040D
233	MONSANTO BLUE C	A0126	16.5 MICRONS	ROUND, SEMI-DULL, WHITE	TYPE L
293	SHAKESPEARE CO. ---•---	A0400	---•---	ROUND, CLEAR, WHITE	WP-89
276	TOLARAM FIBERS INC. ---•---	A0331	---•---	ROUND, SEMI-DULL, WHITE	---•---

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
277	TOLARAM FIBERS INC. ---	A0332	---	ROUND, SEMI-DULL, WHITE	---
278	TOLARAM FIBERS INC. ---	A0333	---	ROUND, CLEAR BRIGHT, WHITE	---
302	WELLMAN, INC. FILLWELL	A0446	---	ROUND, SEMI-DULL, WHITE	TYPE 61G
303	WELLMAN, INC. FILLWELL	K0447	---	ROUND, BRIGHT, GREEN	TYPE 15G
301	WELLMAN, INC. WELLENE	A0445	---	TRILOBAL, BRIGHT, WHITE	TYPE 15BT
298	YKK (USA) INC. NONE	A0424	---	ROUND, SEMI-DULL, WHITE	NONE
299	YKK (USA) INC. NONE	A0425	---	ROUND, SEMI-DULL, WHITE	NONE

POLYESTER/CELLULOSE ACETATES

342	CELANESE LANESE	A0293	23.0 MICRONS	CORE/SHEATH, DULL, WHITE	LANESE-1
343	CELANESE LANESE	A0294	31.2 MICRONS	CORE/SHEATH, DULL, WHITE	LANESE-2

RAYONS

304	AVTEX AVRIL	A0088	10-14 MICRONS	MULTILOBED, ---, WHITE	NO CODE
305	AVTEX AVRIL	A0089	12-16 MICRONS	MULTILOBED, ---, WHITE	NO CODE
306	AVTEX AVRIL	A0092	14-16 MICRONS	BEAN, BRIGHT, WHITE	SN 2656
307	AVTEX AVRIL	A0093	10-15 MICRONS	MULTILOBED, SEMI-DULL, WHITE	SN 2695
308	AVTEX AVRIL	A0095	10-17 MICRONS	MULTILOBED, SEMI-DULL, WHITE	SN 1934
309	AVTEX AVRIL	A0096	25-36 MICRONS	MULTILOBED, BRIGHT, WHITE	SN 2666
310	AVTEX AVRIL	A0097	22-38 MICRONS	MULTILOBED, BRIGHT, WHITE	SN 2659
323	AVTEX DURVIL RAYON	A0310	---	ROUND, BRIGHT, WHITE	SNE 6499
322	AVTEX F-40 HWD RAYON	A0309	13.5 MICRONS	BEAN, BRIGHT, WHITE	SN 2656
325	AVTEX HWM RAYON	A0312	13.5 MICRONS	MULTILOBAL, DULL, WHITE	SN 2695
321	AVTEX REGULAR RAYON	A0308	13.5 MICRONS	ROUND-CRENULATED , BRIGHT, WHITE	SN 2988
324	AVTEX SXD RAYON	A0311	13.5 MICRONS	MULTILOBAL, SEMI-EXTRA DULL, WHITE	SN1934
330	COURTAULDS COLORAY	C0322	13.5 MICRONS	ROUND, BLACK, BLACK	3012
326	COURTAULDS FIBRO	A0317	19 MICRONS	ROUND, BRIGHT, WHITE	1134
327	COURTAULDS FIBRO	A0318	31.1 MICRONS	ROUND, BRIGHT, WHITE	21352
328	COURTAULDS FIBRO	A0320	13.5 MICRONS	ROUND, SEMI-EXTRA DULL, WHITE	18452
329	COURTAULDS FIBRO	A0321	19 MICRONS	ROUND, DULL, WHITE	15432
312	ENKA (BASF) ABSORBIT	A0171	20.0 MICRONS	REGULAR RAYON, DULL, WHITE	8287
311	ENKA (BASF) I. T.	A0168	13.5 MICRONS	REG RAYON-CRENULATED, BRIGHT, WHITE	8644
317	NORTH AMERICAN RAYON BEAU GRIP	F0228	13.5 MICRONS	KIDNEY, DARK, BROWN	9242
319	NORTH AMERICAN RAYON FLOCKING TOW	A0230	13.5 MICRONS	ROUND-SERRATED, DULL, WHITE	9823
320	NORTH AMERICAN RAYON FLOCKING TOW	A0231	25.8 MICRONS	ROUND-SERRATED, BRIGHT, WHITE	9817
313	NORTH AMERICAN RAYON NO TRADE NAME	A0223	15.2 MICRONS	ROUND-SERRATED, EXTRA-DULL, WHITE	3212

<u>No.</u>	<u>Producer/Trade Name</u>	<u>Fiber No.</u>	<u>Diameter</u>	<u>Appearance</u>	<u>Manuf. Code</u>
314	NORTH AMERICAN RAYON NO TRADE NAME	A0225	34.8 MICRONS	FLAT, BRIGHT, WHITE	5815
315	NORTH AMERICAN RAYON NO TRADE NAME	A0226	26.9 MICRONS	ROUND-SERRATED, BRIGHT, WHITE	6412
318	NORTH AMERICAN RAYON T-45	A0229	13.5 MICRONS	ROUND-SERRATED, BRIGHT, WHITE	9813
316	NORTH AMERICAN RAYON TYPE 3000	A0227	13.5 MICRONS	KIDNEY, BRIGHT, WHITE	8707
SARANS					
331	AMETEK NO TRADE NAME	A0152	.71mm W 57 MICRONS TH	FLAT, OPAQUE, WHITE	0514317A
332	AMETEK NO TRADE NAME	A0153	203 MICRONS	ROUND, TRANSPARENT, WHITE	05168317A
SULFAR					
344	PHILLIPS RYTON	A0327	----	----•---, ----•---, WHITE	RYTON
TRIACETATE/POLYESTER					
341	CELANES ARNEL/POLYESTER	A0291	16.0 MICRONS	TRILOBAL, DULL, WHITE	D-161
VINYON					
333	AVTEX NO TRADE NAME	A0098	16-23 MICRONS	ROUND, BRIGHT/SEMI-DULL/DULL, WHITE	HH-2