



HPLC Instrument IQ Version 5.7
Product Release Notes
for Chromeleon 7

Thermo Scientific Chromeleon

Date: September 2022

© 2022 Thermo Fisher Scientific Inc.

1 Overview

HPLC Instrument IQ, version 5.7 is released for Thermo Scientific™ Chromeleon™ Chromatography Data System.

Remark:

The HPLC Instrument IQ Operating Instructions version 5.7 become effective.

All trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries.

2 Changes

2.1 Support of New Thermo Scientific Vanquish Modules [300406]

HPLC Instrument IQ version 5.7 adds support for the following new Thermo Scientific™ Vanquish™ Flex module:

Model Type	Model Number(s)
Vanquish Fraction Collector	VF-F20-A

The instrument method and report files have been updated accordingly.

2.2 Resolved Issue

This section describes the incident that has been resolved with the release of HPLC Instrument IQ version 5.7. The ID column in the table shows the Thermo Fisher Scientific incident IDs.

ID	Description
321163	ISQ detector: Support of new model number "ISQ ICMS Family" introduced with ISQ EC/EM driver update – version 5.0 SP1.