

ation Wizard

You will



you have all
your licens
our compute
e *Install* butt
install it(the

Install

ID: XXRCD-

mNova License Installation Instructions - Notepad

File Edit Format View Help

To update your mNova license files follow these instructions.

mNova will email you a set of one year license codes. Download the license codes.

Open Mestrenova

Choose Help>License Manager

Click Get/Install Licenses

Click Install, you will be prompted to browse for the license files.

Once the licenses are installed, you'll need to close mNova and restart it.

At this point mNova should recognize the file.

If you have any questions, please don't hesitate to contact support.