

MSQ Plus 2.0 SP5 Release Notes

Features

These release notes briefly list changes in the Thermo MSQ Plus™ instrument control application version 2.0 SP5.

Contents

- [Features](#)
- [Minimum requirements](#)
- [Important Information](#)
- [Resolved Issues](#)
- [Known Issues](#)
- [Trademarks](#)

For information on installing the MSQ Plus 2.0 SP5 software, refer to the important information below. For information on configuring and using the MSQ Plus 2.0 SP5 application, refer to the manuals available as PDF files or the Help.

This release incorporates support for Microsoft™ Windows 10 Enterprise (64-bit) operating system.

Minimum Requirements

These are the minimum hardware and software configurations required for MSQ Plus 2.0 SP5 operation.

System	Requirements
Hardware	<ul style="list-style-type: none">• 2.0 GHz processor with 1 GB RAM• Display resolution 1280 x 1024 (XGA)• 75 GB hard drive capacity in NTFS format
Software	<ul style="list-style-type: none">• Microsoft Windows 10 Enterprise (64-bit)• Thermo Foundation™ 3.0• Adobe™ Reader™ 10.1

Important Information

Note:

- It is required to install MSQ Plus driver on either a brand-new PC or fresh image of operating system.
- It is a pre-requisite to have MSQ 2.0 SP4 installed before installing MSQ 2.0 SP5.

Release of MSQ 2.0 SP4 can be found on Knowledge Base (article # [329277](#)).

Follow these instructions to install the MSQ Plus software:

1. Download the [MSQ Plus 2.0 SP5 driver .zip file](#).
2. Navigate to the downloads and unzip the file.
3. To launch the installer, double-click **setup.exe**.
4. Follow the prompts and install the software.

Table 1 lists the defects that were resolved between the MSQ Plus 2.0 SP4 and MSQ Plus 2.0 SP5 applications.

Resolved Issues	Table 1. Resolved issues.		
	Item ID	Software section	Abstract
	87252	Acquisition	Instrument method issues after installing Windows updates
Trademarks	<p>MSQ Plus and Foundation are trademarks, and Thermo Scientific is a registered trademark of Thermo Fisher Scientific Inc. in the United States.</p> <p>Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and other countries. Adobe and Reader are registered trademarks of Adobe Systems Inc. in the United States and other countries.</p> <p>All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.</p>		