

Product description:

Gauge Glass but also referred to as Sight Glass

Part Number 06038


Sight Glass 5/8 X 9

This is a spare part for 5 gallon and 10 gallon classic stills. The model numbers are as follows:

Model A1015-B-61 (also known as model 6772) 240V, 1 PH 5 gallon

A1015-C (also known as model 6773) 120/208V, 3PH 5 gallon

A1016-D (also known as model 6774) 240V, 3PH 10 gallon



Customer issue or question:

It may be difficult to locate information about part number 06038, because of differing descriptions of the part. It may be called the gage glass, glass gauge, glass gage, glass tube, sight glass. Most often it is called sight glass. The Operating Manual and Parts List for the Classic Electrically Heated Stills calls it the Gauge Glass

Service ticket notes

When searching for information, please recognize that any of these descriptions are correct for part number 06038